



Cholesterol FS*

In-vitro-Diagnostic for veterinary use only

Diagnostic reagent for quantitative in vitro determination of cholesterol in serum or plasma on DiaSys respons[®]910 VET

Order Information


Cat. No. 1 1300 99 11 923
4 containers for 200 tests each

Method

"CHOD-PAP": enzymatic photometric test

Principle

Determination of cholesterol after enzymatic hydrolysis and oxidation. The colorimetric indicator is quinoneimine which is generated from 4-aminoantipyrine and phenol by hydrogen peroxide under the catalytic action of peroxidase (Trinder's reaction).



Reagent

Components and Concentrations

| | | |
|----------------------|--------|------------|
| Good's buffer | pH 6.7 | 50 mmol/L |
| Phenol | | 5 mmol/L |
| 4-Aminoantipyrine | | 0.3 mmol/L |
| Cholesterol esterase | (CHE) | ≥ 200 U/L |
| Cholesterol oxidase | (CHO) | ≥ 50 U/L |
| Peroxidase | (POD) | ≥ 3 kU/L |

Storage Instructions and Reagent Stability

The reagent is stable up to the end of the indicated month of expiry, if stored at 35.6 – 46.4°F, protected from light and contamination is avoided. DiaSys response containers provide protection from light. Do not freeze the reagent!

Warnings and Precautions

1. The reagent contains sodium azide (0.95 g/L) as preservative. Do not swallow! Avoid contact with skin and mucous membranes.
2. In very rare cases, samples of animals with gammopathy might give falsified results.
3. N-acetylcysteine (NAC), acetaminophen and metamizole medication leads to falsely low results.
4. Please refer to the safety data sheets and take the necessary precautions for the use of laboratory reagents. For diagnostic purposes, the results should always be assessed with the animal's medical history, clinical examinations and other findings.
5. For professional use only!

Waste Management

Please refer to local legal requirements.

Reagent Preparation

The reagent is ready to use. The bottles are placed directly onto the reagent rotor.

Specimen

Serum, heparin plasma or EDTA plasma

Stability :
2 days 39.2°F to 46.4°F

Discard contaminated specimens.

Calibrators and Controls

For the calibration the DiaSys TruCal U calibrator is recommended. The assigned values of the calibrator have been made traceable to the reference method gas chromatography-isotope dilution mass spectrometry (GC-IDMS). For internal quality control DiaSys TruLab N and P or TruLab L controls should be assayed. Each laboratory should establish corrective actions in case of deviations in control recovery.

| | Cat. No. | Kit size |
|------------------|------------------|-----------|
| TruCal U | 5 9100 99 11 063 | 20 x 3 mL |
| TruLab N | 5 9000 99 11 062 | 20 x 5 mL |
| TruLab P | 5 9050 99 11 062 | 20 x 5 mL |
| TruLab L Level 1 | 5 9020 99 11 065 | 3 x 3 mL |
| TruLab L Level 2 | 5 9030 99 11 065 | 3 x 3 mL |

Performance Characteristics

The performance characteristics were evaluated with human samples and might differ from results obtained with various animal specimen.

| | |
|---|---------------------|
| Measuring range up to 750 mg/dL cholesterol (in case of higher concentrations re-measure samples after manual dilution with NaCl solution (9 g/L) or use rerun function). | |
| Limit of detection** | 1 mg/dL cholesterol |
| On-board stability | 8 weeks |
| Calibration stability | 4 weeks |

| Interfering substance | Interferences < 10% | Cholesterol [mg/dL] |
|-------------------------|---------------------|---------------------|
| Ascorbate | up to 6 mg/dL | 222 |
| Hemoglobin | up to 230 mg/dL | 152 |
| | up to 230 mg/dL | 223 |
| Bilirubin, conjugated | up to 15 mg/dL | 147 |
| | up to 25 mg/dL | 236 |
| Bilirubin, unconjugated | up to 21 mg/dL | 149 |
| | up to 23 mg/dL | 237 |
| Lipemia (triglycerides) | up to 2200 mg/dL | 136 |
| | up to 2200 mg/dL | 234 |





For further information on interfering substances refer to Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th. ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press, 2000.

** according to NCCLS document EP17-A, vol. 24, no. 34

Conversion Factor

Cholesterol [mg/dL] x 0.02586 = Cholesterol [mmol/L]

Reference Range

|  |  |  |  | Unit |
|---|---|---|---|-------|
| DOG | CAT | HORSE | CATTLE | |
| 124 – 327 | 63 – 258 | 67 – 136 * | 132 – 332 | mg/dL |

Source:

Reference ranges have been validated by DiaSys USA according to National Reference Laboratory standards.

* Estimated: Based on preliminary results and findings in the literature.

Each laboratory should check if the reference ranges are transferable to its own animal population and determine own reference ranges if necessary.

Manufacturer

DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim Germany