



# CK-NAC FS\* IFCC

In-vitro-Diagnostic for veterinary use only

Diagnostic reagent for quantitative in vitro determination of creatine kinase (CK) in serum or plasma on DiaSys respons<sup>®</sup>910 VET


## Order Information

Cat. No. 1 1601 99 11 921  
4 twin containers for 120 tests each

## Method

Optimized UV-test according to IFCC (International Federation of Clinical Chemistry and Laboratory Medicine) and DGKC (German Society of Clinical Chemistry)

## Principle



## Reagents

### Components and Concentrations

<b>R1:</b>	Imidazole	pH 6.0	60 mmol/L
	Glucose		27 mmol/L
	N-Acetylcysteine	(NAC)	27 mmol/L
	Magnesium acetate		14 mmol/L
	EDTA-Na <sub>2</sub>		2 mmol/L
	NADP		2.7 mmol/L
	Hexokinase	(HK)	≥ 5 kU/L
<b>R2:</b>	Imidazole	pH 9.0	160 mmol/L
	ADP		11 mmol/L
	AMP		28 mmol/L
	Diadenosine pentaphosphate		55 μmol/L
	Glucose-6-phosphate dehydrogenase (G6P-DH)		≥ 14 kU/L
	EDTA-Na <sub>2</sub>		2 mmol/L
	Creatine phosphate		160 mmol/L

### Storage Instructions and Reagent Stability

The reagents are stable up to the end of the indicated month of expiry, if stored at 35.6 – 46.4°F, protected from light and contamination is avoided. DiaSys respons containers provide protection from light. Do not freeze the reagents!

### Warnings and Precautions

1. Reagent 1: Danger. H360D May damage the unborn child. P201 Obtain special instructions before use. P280 Wear protective gloves/protective clothing/eye protection/face protection. P308+P313 If exposed or concerned: Get medical advice/attention.
2. Reagent 2: Danger. H315 Causes skin irritation. H319 Causes serious eye irritation. H360D May damage the unborn child. P201 Obtain special instructions before use. P280 Wear protective gloves/protective clothing/eye protection/face protection. P302+P352 If on skin: Wash with plenty of water/soap. P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P308+P313 If exposed or concerned: Get medical advice/attention.
3. The reagents contain sodium azide (0.95 g/L) as preservative. Do not swallow! Avoid contact with skin and mucous membranes.
4. Reagent 2 contains animal material. Handle the product as potentially infectious according to universal precautions and good laboratory practice.
5. In very rare cases, samples of animals with gammopathy might give falsified results.
6. Please refer to the safety data sheets and take the necessary precautions for the use of laboratory reagents. For diagnostic purposes, the results should always be assessed with the animal's medical history, clinical examinations and other findings.
7. For professional use only!

### Waste Management

Please refer to local legal requirements.

### Reagent Preparation

The reagents are ready to use. The bottles are placed directly onto the reagent rotor.

### Specimen

Serum, heparin plasma or EDTA plasma

Stability :  
2 days 39.2°F to 46.4°F

Discard contaminated specimens.

## Calibrators and Controls

For calibration, DiaSys TruCal U calibrator is recommended. This method has been standardized against the original IFCC formulation. For internal quality control DiaSys TruLab N and P controls should be assayed. Each laboratory should establish corrective action in case of deviations in control recovery.

	Cat. No.	Kit size
TruCal U	5 9100 99 11 063	20 x 3 mL
TruLab N	5 9000 99 11 062	20 x 5 mL
TruLab P	5 9050 99 11 062	20 x 5 mL

## Performance Characteristics

The performance characteristics were evaluated with human samples and might differ from results obtained with various animal specimen.

Measuring range up to 1100 U/L CK (in case of higher activities re-measure samples after manual dilution with NaCl solution (9 g/L) or use rerun function).	
Limit of detection**	3 U/L CK
On-board stability	6 weeks
Calibration stability	3 weeks

Interfering substance	Interferences < 10%	CK-NAC [U/L]
Ascorbate	up to 30 mg/dL	99.0
Hemoglobin	up to 100 mg/dL	143
	up to 100 mg/dL	197
Bilirubin, conjugated	up to 60 mg/dL	92.0
	up to 60 mg/dL	175
Bilirubin, unconjugated	up to 70 mg/dL	96.7
	up to 70 mg/dL	307
Lipemia (triglycerides)	up to 1000 mg/dL	90.5
	up to 2000 mg/dL	158





For further information on interfering substances refer to Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th. ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press, 2000.

\*\* according to NCCLS document EP17-A, vol. 24, no. 34

## Conversion Factor

CK-NAC [U/L] x 0.0167 = CK-NAC [μkat/L]

## Reference Range

				Unit
DOG	CAT	HORSE	CATTLE	
57 – 308	79 – 778	185 – 470	130 – 575	U/L

Source:

Reference ranges have been validated by DiaSys USA according to National Reference Laboratory standards.

Each laboratory should check if the reference ranges are transferable to its own animal population and determine own reference ranges if necessary. For diagnostic purposes CK values should always be assessed in conjunction with the anamnesis, the clinical examination and other findings.

## Manufacturer

DiaSys Diagnostic Systems GmbH  
Alte Strasse 9 65558 Holzheim Germany