



TruCal U

In-vitro-Diagnostic for veterinary use only

Multi-calibrator for use in tests for quantitative in vitro determination of various analytes on DiaSys respons[®]910 VET

Order Information

Cat. No.	Kit size
5 9100 99 11 063	20 x 3 mL
5 9100 99 11 064	6 x 3 mL

Description

TruCal U is a lyophilized calibrator based on human blood material (serum) and contains chemical additives and biological material of specified origin.

Analyte	Origin
Acid phosphatase total	Human prostate/potato
Albumin	Bovine plasma
Aldolase	Rabbit muscle
Alkaline phosphatase	Placenta (human, recombinant)
Alanine-Aminotransferase	Porcine (heart)
α-Amylase	Porcine (pancreas)
Pancreatic amylase	Porcine (pancreas)
Aspartate Aminotransferase	Human, recombinant
Cholesterol	Bovine plasma
Cholinesterase	Human serum
Creatin kinase	Rabbit muscle
γ-Glutamyltransferase	Human, recombinant
GLDH	Bacterium, recombinant
Lactatedehydrogenase	Porcine heart
Lipase	Pancreas (human recombinant)
Triglycerides	Chicken egg yolk

The concentration of the biological material does not exceed the maximum, lot specific target value concentration of the analyte.

Storage and Stability

Unopened: Until the end of the indicated month of expiry, if stored at 35.6 – 46.4°F.

Once reconstituted, TruCal U calibrator can be used within the period reported in the table below if stored tightly closed at the indicated temperature and bacterial contamination is avoided.

Stability of the components after reconstitution:

	-4°F *	+39.2°F	+77°F
Bilirubin (kept in the dark)	14 days	8 hours	3 hours
Other analytes	30 days	2 days	8 hours

* freeze only once!

Criterion for the stability stated is a recovery within ±5% of the initial value.

Warnings and Precautions

- Each individual blood donation used for production of TruCal U was found to be non-reactive when tested with FDA-approved methods or CE-accepted methods for HBsAG, anti-HIV 1+2 and anti-HCV. As there is no possibility to exclude definitely that products derived from human blood or animal components transmit infectious agents, it is recommended to handle TruCal U calibrator with the same precautions used for specimens.
- TruCal U contains biological material of specified origin. The calibrators should be handled as potentially infectious and with the same precautions used for specimens.
- Please refer to the safety data sheets and take the necessary precautions for the use of calibrators and controls.
- For professional use only!

Preparation

The lyophilizate is sealed under vacuum; therefore the vial should be opened very carefully to avoid loss of dried material. For reconstitution add exactly 3 mL of distilled water. Close the vial carefully and allow the calibrator to stand for 30 minutes swirling occasionally. Avoid foaming! Do not shake!

Defrost reconstituted, frozen aliquots of TruCal U protected from light at room temperature (64.4 – 77°F). To homogenize after complete defrosting, slightly swivel aliquots and immediately afterwards use them for calibration in the same way as the freshly reconstituted TruCal U.

Procedure

Please refer to the reagent package inserts for instructions for use.

Calibrator Values

The concentrations of the calibrator analytes are lot-specific and given in the value sheet of the corresponding lot. The values were determined using the method and using the reagents specified by the given catalogue number. Determinations were performed under standardized protocols using DiaSys reagents and TruCal U master calibrator or reference materials. The value sheet contains information on traceability.

Waste Management

Please refer to local legal requirements.

Manufacturer

DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim Germany