



Cleaning solutions for use on DiaSys respons[®]910 VET

In-vitro-Diagnostic for veterinary use only

Cleaner A

Acidic cleaning solution

Order Information

Cat. No. 1 8610 99 11 923

4 containers with 60 mL each

Storage Instructions and Stability

Cleaner A is stable up to the end of the indicated month of expiry if stored at 35.6 – 77°F and contamination is avoided. Do not freeze!

Warnings and Precautions

1. Warning. H290 May be corrosive to metals. P262 Do not get in eyes, on skin, or on clothing. P390 Absorb spillage to prevent material damage. P406 Store in corrosive resistant plastic container/in container with a resistant inner liner.
2. Please refer to the safety data sheet and take the necessary precautions for the use of laboratory reagents.
3. For professional use only!

Cleaner B

Alkaline cleaning solution

Order Information

Cat. No. 1 8650 99 11 923

4 containers with 60 mL each

Storage Instructions and Stability

Cleaner B is stable up to the end of the indicated month of expiry if stored at 35.6 – 77°F and contamination is avoided. Do not freeze!

Warnings and Precautions

1. Warning. H290 May be corrosive to metals. P234 Keep only in original container. P390 Absorb spillage to prevent material damage.
2. Please refer to the safety data sheet and take the necessary precautions for the use of laboratory reagents.
3. For professional use only!

Manufacturer

DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim Germany