



Cleaning solutions for use on DiaSys respons[®]910 VET

In-vitro-Diagnostic for veterinary use only

Cleaner A

Acidic cleaning solution

Order Information

Cat. No. 1 8610 99 11 923

4 containers with 60 mL each

Storage Instructions and Stability

Cleaner A is stable up to the end of the indicated month of expiry if stored at 35.6 – 77°F and contamination is avoided. Do not freeze!

Warnings and Precautions

- Warning. H290 May be corrosive to metals. H315 Causes skin irritation. H319 Causes serious eye irritation. P234 Keep only in original container. P280 Wear protective gloves/protective clothing/eye protection/face protection. P302+P352 If on skin: Wash with plenty of water/soap. P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P337+P313 If eye irritation persists: Get medical advice/attention. P390 Absorb spillage to prevent material damage.
- Please refer to the safety data sheet and take the necessary precautions for the use of laboratory reagents.
- For professional use only!

Cleaner B

Alkaline cleaning solution

Order Information

Cat. No. 1 8650 99 11 923

4 containers with 60 mL each

Storage Instructions and Stability

Cleaner B is stable up to the end of the indicated month of expiry if stored at 35.6 – 77°F and contamination is avoided. Do not freeze!

Warnings and Precautions

- Warning. H290 May be corrosive to metals. P234 Keep only in original container. P390 Absorb spillage to prevent material damage.
- Please refer to the safety data sheet and take the necessary precautions for the use of laboratory reagents.
- For professional use only!

Manufacturer

DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim Germany