



# SAFETY DATA SHEET

in accordance with 29 CFR 1910.1200 and ANSI standard Z400.1-2010

## NEFA FS Reagent R2

Material number 1 5781 R2

Revision date: 11/19/2018

Version: 12

Language: en-US

Date of print: 1/30/2019

Page: 1 of 6

## 1. Product and company identification

### Product identifier

Trade name: NEFA FS Reagent R2  
As part of the kits: 1 5781 XX XX XXX  
(The positions X code different packages.)

### Relevant identified uses of the substance or mixture and uses advised against

General use: Reagent for in-vitro diagnostic use  
For professional use only

### Details of the supplier of the safety data sheet

Company name: DiaSys Diagnostic Systems GmbH  
Street/POB-No.: Alte Strasse 9  
Postal Code, city: 65558 Holzheim  
WWW: <http://www.diasys.de>  
E-mail: [mail@diasys.de](mailto:mail@diasys.de)  
Telephone: +49 (0) 6432-9146-0  
Telefax: +49 (0) 6432-9146-32  
Dept. responsible for information:  
Corporate headquarters, Telephone: +49 (0) 6432-9146-0, Email: [mail@diasys.de](mailto:mail@diasys.de)

### Emergency phone number

Infraserv, Telephone: +49 (0) 69-305-6418

## 2. Hazards identification

### Emergency overview

Appearance: Physical state at 68 °F and 101.3 kPa: liquid  
Color: yellow up to brownish, clear  
Odor: odorless  
Classification: Eye Irritation - Category 2A;

Hazard symbols:



Signal word: **Warning**  
Hazard statements: Causes serious eye irritation.  
Precautionary statements:

Wear protective gloves/protective clothing/eye protection/face protection.  
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
If eye irritation persists: Get medical advice/attention.

### Regulatory status

This material is considered hazardous by the U.S. OSHA Hazard Communication Standard (29 CFR 1910.1200) and SIMDUT in Canada.

### Hazards not otherwise classified

see section 11: Toxicological information



# SAFETY DATA SHEET

in accordance with 29 CFR 1910.1200 and ANSI standard Z400.1-2010

## NEFA FS Reagent R2

Material number 1 5781 R2

Revision date: 11/19/2018

Version: 12

Language: en-US

Date of print: 1/30/2019

Page: 2 of 6

### 3. Composition / Information on ingredients

Chemical characterization: Aqueous solution of anorganic salts and organic compounds.

Relevant ingredients:

CAS No.	Designation	Content	Classification
CAS 9002-93-1	Octylphenoxy polyethoxyethanol	< 2.5 %	Acute Toxicity - oral - Category 4. Eye Damage - Category 1. Aquatic toxicity - chronic - Category 2.
CAS 10191-18-1	2-[bis(2-Hydroxyethyl)amino] ethanesulphonic acid	< 2 %	Skin Irritation - Category 2. Eye Irritation - Category 2A. Specific Target Organ Toxicity (Single Exposure) - Category 3.

### 4. First aid measures

General information:	First aider: Pay attention to self-protection! If medical advice is needed, have product container or label at hand.
In case of inhalation:	Provide fresh air. Seek medical aid in case of troubles.
Following skin contact:	Change contaminated clothing. Remove residues with water. Seek medical aid in case of troubles.
After eye contact:	Immediately flush eyes with plenty of flowing water for 10 to 15 minutes holding eyelids apart. Remove contact lenses, if present and easy to do. Continue rinsing. Subsequently seek the immediate attention of an ophthalmologist.
After swallowing:	Rinse mouth immediately and drink plenty of water. Never give anything by mouth to an unconscious person. Seek medical attention.

#### Most important symptoms/effects, acute and delayed

Causes serious eye irritation.

#### Information to physician

Treat symptomatically.

### 5. Fire fighting measures

Flash point/flash point range:

not combustible

Auto-ignition temperature: No data available

Suitable extinguishing media:

Product is non-combustible. Extinguishing materials should therefore be selected according to surroundings.

#### Specific hazards arising from the chemical

Fires in the immediate vicinity may cause the development of dangerous vapors. In case of fire may be liberated: Nitrogen oxides (NO<sub>x</sub>), sulphur oxides, Carbon monoxide and carbon dioxide.

Protective equipment and precautions for firefighters:

Wear self-contained breathing apparatus.

Additional information: Do not allow fire water to penetrate into surface or ground water.



# SAFETY DATA SHEET

in accordance with 29 CFR 1910.1200 and ANSI standard Z400.1-2010

## NEFA FS Reagent R2

Material number 1 5781 R2

Revision date: 11/19/2018

Version: 12

Language: en-US

Date of print: 1/30/2019

Page: 3 of 6

### 6. Accidental release measures

- Personal precautions: Avoid contact with skin and eyes. Wear appropriate protective equipment. Provide adequate ventilation. Keep unprotected people away.
- Environmental precautions: Do not allow to penetrate into soil, waterbodies or drains.
- Methods for clean-up: Soak up with absorbent materials such as sand, siliceus earth, acid- or universal binder. Store in special closed containers and dispose of according to ordinance. Wash spill area with plenty of water.

### 7. Handling and storage

#### Handling

Advices on safe handling: Avoid contact with skin and eyes. Wear appropriate protective equipment. Keep all containers, equipment and working place clean. Provide adequate ventilation, and local exhaust as needed. Do not eat, drink or smoke when using this product. Wash hands thoroughly after handling. Have eye wash bottle or eye rinse ready at work place.

#### Storage

Requirements for storerooms and containers:

Keep containers tightly closed and at a temperature between 35.6 °F and 46.4 °F. Do not freeze. Protect from light. Keep sterile.

Hints on joint storage: Keep away from food, drink and animal feedingstuffs.

### 8. Exposure controls / personal protection

#### Engineering controls

Provide good ventilation and/or an exhaust system in the work area.  
See also information in chapter 7, section storage.

#### Personal protection equipment (PPE)

- Eye/face protection: Tightly sealed goggles according to OSHA Standard - 29 CFR: 1910.133 or ANSI Z87.1-2010.
- Skin protection: Wear suitable protective clothing.  
Protective gloves according to OSHA Standard - 29 CFR: 1910.138.  
Glove material: Butyl caoutchouc (butyl rubber) - Breakthrough time: >480 min.  
Observe glove manufacturer's instructions concerning penetrability and breakthrough time.
- Respiratory protection: Use a breathing protection against vapors/aerosol.  
Use combination filter type A/P according to OSHA Standard - 29 CFR: 1910.134 or ANSI Z88.2.
- General hygiene considerations: Avoid contact with skin and eyes. Change contaminated clothing. Wash hands before breaks and after work. Do not eat, drink or smoke when using this product. Have eye wash bottle or eye rinse ready at work place.



# SAFETY DATA SHEET

in accordance with 29 CFR 1910.1200 and ANSI standard Z400.1-2010

## NEFA FS Reagent R2

Material number 1 5781 R2

Revision date: 11/19/2018

Version: 12

Language: en-US

Date of print: 1/30/2019

Page: 4 of 6

## 9. Physical and chemical properties

### Information on basic physical and chemical properties

Appearance:	Physical state at 68 °F and 101.3 kPa: liquid Color: yellow up to brownish, clear
Odor:	odorless
Odor threshold:	No data available
pH value:	at 77 °F: 7.0
Melting point/freezing point:	No data available
Initial boiling point and boiling range:	No data available
Flash point/flash point range:	not combustible
Evaporation rate:	No data available
Flammability:	No data available
Explosion limits:	No data available
Vapor pressure:	No data available
Vapor density:	No data available
Density:	at 68 °F: 1.0075 g/mL
Water solubility:	completely miscible
Partition coefficient: n-octanol/water:	No data available
Auto-ignition temperature:	No data available
Thermal decomposition:	No data available
Additional information:	No data available

## 10. Stability and reactivity

Reactivity:	No data available
Chemical stability:	Stable under recommended storage conditions.
Possibility of hazardous reactions	No hazardous reactions known.
Conditions to avoid:	Protect against heat /sun rays.
Incompatible materials:	No data available
Hazardous decomposition products:	No hazardous decomposition products when regulations for storage and handling are observed.
Thermal decomposition:	No data available



# SAFETY DATA SHEET

in accordance with 29 CFR 1910.1200 and ANSI standard Z400.1-2010

## NEFA FS Reagent R2

Material number 1 5781 R2

Revision date: 11/19/2018

Version: 12

Language: en-US

Date of print: 1/30/2019

Page: 5 of 6

## 11. Toxicological information

### Toxicological tests

Toxicological effects: The statements are derived from the properties of the single components. No toxicological data is available for the product as such.

Acute toxicity (oral): Lack of data.

Acute toxicity (dermal): Lack of data.

Acute toxicity (inhalative): Lack of data.

Skin corrosion/irritation: Lack of data.

Serious eye damage/irritation: Eye Irritation - Category 2A = Causes serious eye irritation.

Sensitisation to the respiratory tract: Lack of data.

Skin sensitisation: Lack of data.

Germ cell mutagenicity/Genotoxicity: Lack of data.

Carcinogenicity: Lack of data.

Reproductive toxicity: Lack of data.

Effects on or via lactation: Lack of data.

Specific target organ toxicity (single exposure): Lack of data.

Specific target organ toxicity (repeated exposure): Lack of data.

Aspiration hazard: Lack of data.

## 12. Ecological information

### Ecotoxicity

Further details: No data available

### Mobility in soil

No data available

### Persistence and degradability

Further details: No data available

### Additional ecological information

General information: Do not allow to enter into ground-water, surface water or drains.

## 13. Disposal considerations

### Product

Recommendation: Special waste. Dispose of waste according to applicable legislation.

### Contaminated packaging

Recommendation: Dispose of waste according to applicable legislation.  
Non-contaminated packages may be recycled.



# SAFETY DATA SHEET

in accordance with 29 CFR 1910.1200 and ANSI standard Z400.1-2010

## NEFA FS Reagent R2

Material number 1 5781 R2

Revision date: 11/19/2018

Version: 12

Language: en-US

Date of print: 1/30/2019

Page: 6 of 6

### 14. Transport information

#### USA: Department of Transportation (DOT)

Proper shipping name: Not restricted

#### Sea transport (IMDG)

Proper shipping name: Not restricted

Marine pollutant: no

#### Air transport (IATA)

Proper shipping name: Not restricted

#### Further information

No dangerous good in sense of these transport regulations.

### 15. Regulatory information

#### National regulations - U.S. Federal Regulations

Octylphenoxypolyethoxyethanol: TSCA Inventory: listed; EPA flags XU

TSCA HPVC: not listed

2-[bis(2-Hydroxyethyl)amino]ethanesulphonic acid: TSCA Inventory: listed

TSCA HPVC: not listed

#### National regulations - Great Britain

Hazchem-Code: -

### 16. Other information

Text for labeling: Contains < 2.5 % Octylphenoxypolyethoxyethanol, < 2 %  
2-[bis(2-Hydroxyethyl)amino]ethanesulphonic acid. Safety data sheet available on request.

Hazard rating systems: NFPA Hazard Rating:



Health: 1 (Slight)  
Fire: 0 (Minimal)  
Reactivity: 0 (Minimal)

HMIS Version III Rating:

Health: 1 (Slight)  
Flammability: 0 (Minimal)  
Physical Hazard: 0 (Minimal)

Personal Protection: X = Consult your supervisor

HEALTH	1
FLAMMABILITY	0
PHYSICAL HAZARD	0
	X

Reason of change: Changes in section 2: Labeling

Date of first version: 3/13/2008

#### Department issuing data sheet

Contact person: see section 1: Dept. responsible for information

The information in this data sheet has been established to our best knowledge and was up-to-date at time of revision. It does not represent a guarantee for the properties of the product described in terms of the legal warranty regulations.