



# SAFETY DATA SHEET

in accordance with 29 CFR 1910.1200 and ANSI standard Z400.1-2010

Revision date: 8/28/2018

Version: 1

Language: en-US

Date of print: 1/30/2019

## $\beta$ -Hydroxybutyrate 21 Tri-Level Controls, $\beta$ -Hydroxybutyrate 21 Bi-Level Controls

Material number 5 9570/5 9540

Page: 1 of 6

### 1. Product and company identification

#### Product identifier

Trade name:  $\beta$ -Hydroxybutyrate 21 Tri-Level Controls,  $\beta$ -Hydroxybutyrate 21 Bi-Level Controls  
As part of the kits: 5 9570 XX XX XXX (Tri-Level); 5 9540 XX XX XXX (Bi-Level)  
(The positions X code different packages.)

#### Relevant identified uses of the substance or mixture and uses advised against

General use: Reagent for in-vitro diagnostic use

#### Details of the supplier of the safety data sheet

Company name: DiaSys Diagnostic Systems GmbH

Street/POB-No.: Alte Strasse 9

Postal Code, city: 65558 Holzheim

WWW: <http://www.diasys.de>

E-mail: [mail@diasys.de](mailto:mail@diasys.de)

Telephone: +49 (0) 6432-9146-0

Telefax: +49 (0) 6432-9146-32

Dept. responsible for information:

Corporate headquarters, Telephone: +49 (0) 6432-9146-0, Email: [mail@diasys.de](mailto:mail@diasys.de)

#### Emergency phone number

Infraserv, Telephone: +49 (0) 69-305-6418

### 2. Hazards identification

#### Emergency overview

Appearance: Physical state at 68 °F and 101.3 kPa: liquid

Color: colorless

Odor: odorless

Classification: This material is classified as not hazardous.

#### Regulatory status

This material is not considered hazardous by the U.S. OSHA Hazard Communication Standard (29 CFR 1910.1200) and SIMDUT in Canada.

#### Hazards not otherwise classified

see section 11: Toxicological information

### 3. Composition / Information on ingredients

Chemical characterization: Aqueous solution of anorganic salts and organic compounds.

Additional information: Preparation does not contain dangerous substances above limits that need to be mentioned in this section according to applicable legislation.



# SAFETY DATA SHEET

in accordance with 29 CFR 1910.1200 and ANSI standard Z400.1-2010

Revision date: 8/28/2018

Version: 1

Language: en-US

Date of print: 1/30/2019

## $\beta$ -Hydroxybutyrate 21 Tri-Level Controls, $\beta$ -Hydroxybutyrate 21 Bi-Level Controls

Material number 5 9570/5 9540

Page: 2 of 6

### 4. First aid measures

- In case of inhalation: Move victim to fresh air. In case of respiratory difficulties seek medical attention.
- Following skin contact: Change contaminated clothing. Remove residues with water. In case of skin reactions, consult a physician.
- After eye contact: Immediately flush eyes with plenty of flowing water for 10 to 15 minutes holding eyelids apart. Remove contact lenses, if present and easy to do. Continue rinsing. In case of troubles or persistent symptoms, consult an ophthalmologist.
- After swallowing: Rinse mouth and drink large quantities of water. If you feel unwell, seek medical advice. Never give anything by mouth to an unconscious person.

#### Most important symptoms/effects, acute and delayed

No data available

#### Information to physician

Treat symptomatically.

### 5. Fire fighting measures

- Flash point/flash point range: No data available
- Auto-ignition temperature: No data available
- Suitable extinguishing media: Product is non-combustible. Extinguishing materials should therefore be selected according to surroundings.

#### Specific hazards arising from the chemical

Fires in the immediate vicinity may cause the development of dangerous vapors. sodium oxides, carbon monoxide and carbon dioxide.

- Protective equipment and precautions for firefighters: Wear self-contained breathing apparatus.
- Additional information: Do not allow fire water to penetrate into surface or ground water.

### 6. Accidental release measures

- Personal precautions: Avoid contact with skin and eyes. Wear appropriate protective equipment. Provide adequate ventilation.
- Environmental precautions: Do not allow to penetrate into soil, waterbodies or drains.
- Methods for clean-up: Soak up with absorbent materials such as sand, siliceus earth, acid- or universal binder. Store in special closed containers and dispose of according to ordinance. Wash spill area with plenty of water.



# SAFETY DATA SHEET

in accordance with 29 CFR 1910.1200 and ANSI standard Z400.1-2010

Revision date: 8/28/2018

Version: 1

Language: en-US

Date of print: 1/30/2019

## $\beta$ -Hydroxybutyrate 21 Tri-Level Controls, $\beta$ -Hydroxybutyrate 21 Bi-Level Controls

Material number 5 9570/5 9540

Page: 3 of 6

### 7. Handling and storage

#### Handling

Advices on safe handling: For use only by qualified technical staff. Provide adequate ventilation, and local exhaust as needed. Wear appropriate protective equipment. Avoid contact with skin, eyes, and clothing. After worktime and during work intervals the affected skin areas must be thoroughly cleaned. Take off contaminated clothing and wash it before reuse. Provide a conveniently located eye rinse station.

#### Storage

Requirements for storerooms and containers:

Keep containers tightly closed and at a temperature between 35.6 °F and 46.4 °F. Do not freeze. Protect from light. Keep sterile.

### 8. Exposure controls / personal protection

#### Engineering controls

Provide adequate ventilation, and local exhaust as needed.

See also information in chapter 7, section storage.

#### Personal protection equipment (PPE)

Eye/face protection: Tightly sealed goggles according to OSHA Standard - 29 CFR: 1910.133 or ANSI Z87.1-2010.

Skin protection: Wear suitable protective clothing.

Protective gloves according to OSHA Standard - 29 CFR: 1910.138.

Glove material: Nitrile rubber - Breakthrough time: >480 min.

Observe glove manufacturer's instructions concerning penetrability and breakthrough time.

Respiratory protection: Provide adequate ventilation.

General hygiene considerations:

Avoid contact with skin, eyes, and clothing. Take off contaminated clothing and wash it before reuse. Wash hands before breaks and after work. Provide a conveniently located eye rinse station.

### 9. Physical and chemical properties

#### Information on basic physical and chemical properties

Appearance: Physical state at 68 °F and 101.3 kPa: liquid

Color: colorless

Odor: odorless

Odor threshold: No data available

pH value: at 77 °F: 5.6

Melting point/freezing point: No data available

Initial boiling point and boiling range: No data available

Flash point/flash point range: No data available

Evaporation rate: No data available

Flammability: No data available



# SAFETY DATA SHEET

in accordance with 29 CFR 1910.1200 and ANSI standard Z400.1-2010

Revision date: 8/28/2018

Version: 1

Language: en-US

Date of print: 1/30/2019

## $\beta$ -Hydroxybutyrate 21 Tri-Level Controls, $\beta$ -Hydroxybutyrate 21 Bi-Level Controls

Material number 5 9570/5 9540

Page: 4 of 6

Explosion limits:	No data available
Vapor pressure:	No data available
Vapor density:	No data available
Density:	at 68 °F: 1.02 g/mL
Water solubility:	completely miscible
Partition coefficient: n-octanol/water:	No data available
Auto-ignition temperature:	No data available
Thermal decomposition:	No data available
Additional information:	No data available

### 10. Stability and reactivity

Reactivity:	No data available
Chemical stability:	Stable under recommended storage conditions.
Possibility of hazardous reactions	No hazardous reactions known.
Conditions to avoid:	Protect against heat /sun rays. Protect from frost.
Incompatible materials:	No data available
Hazardous decomposition products:	No decomposition when used properly.
Thermal decomposition:	No data available

### 11. Toxicological information

#### Toxicological tests

Toxicological effects:	Acute toxicity (oral): Lack of data.
	Acute toxicity (dermal): Lack of data.
	Acute toxicity (inhalative): Lack of data.
	Skin corrosion/irritation: Lack of data.
	Serious eye damage/irritation: Lack of data.
	Sensitisation to the respiratory tract: Lack of data.
	Skin sensitisation: Lack of data.
	Germ cell mutagenicity/Genotoxicity: Lack of data.
	Carcinogenicity: Lack of data.
	Reproductive toxicity: Lack of data.
	Effects on or via lactation: Lack of data.
	Specific target organ toxicity (single exposure): Lack of data.
	Specific target organ toxicity (repeated exposure): Lack of data.
	Aspiration hazard: Lack of data.



# SAFETY DATA SHEET

in accordance with 29 CFR 1910.1200 and ANSI standard Z400.1-2010

Revision date: 8/28/2018

Version: 1

Language: en-US

Date of print: 1/30/2019

## $\beta$ -Hydroxybutyrate 21 Tri-Level Controls, $\beta$ -Hydroxybutyrate 21 Bi-Level Controls

Material number 5 9570/5 9540

Page: 5 of 6

### 12. Ecological information

#### Ecotoxicity

Further details: No data available

#### Mobility in soil

No data available

#### Persistence and degradability

Further details: No data available

#### Additional ecological information

General information: Do not allow to enter into ground-water, surface water or drains.

### 13. Disposal considerations

#### Product

Recommendation: Special waste. Dispose of waste according to applicable legislation.

#### Contaminated packaging

Recommendation: Dispose of waste according to applicable legislation. Non-contaminated packages may be recycled.

### 14. Transport information

#### USA: Department of Transportation (DOT)

Proper shipping name: Not restricted

#### Sea transport (IMDG)

Proper shipping name: Not restricted

Marine pollutant: no

#### Air transport (IATA)

Proper shipping name: Not restricted

#### Further information

No dangerous good in sense of these transport regulations.

### 15. Regulatory information

#### National regulations - U.S. Federal Regulations

Trisodium citrate - 2-hydrate: TSCA Inventory: not listed  
TSCA listed: Trisodium citrate, anhydrous CAS No. 68-04-2

Albumins, blood serum: TSCA Inventory: listed; UVCB; EPA flags XU  
TSCA HPVC: not listed



# SAFETY DATA SHEET

in accordance with 29 CFR 1910.1200 and ANSI standard Z400.1-2010

Revision date: 8/28/2018

Version: 1

Language: en-US

Date of print: 1/30/2019

## $\beta$ -Hydroxybutyrate 21 Tri-Level Controls, $\beta$ -Hydroxybutyrate 21 Bi-Level Controls

Material number 5 9570/5 9540

Page: 6 of 6

### National regulations - Great Britain

Hazchem-Code: -

## 16. Other information

Hazard rating systems:



NFPA Hazard Rating:

Health: 1 (Slight)

Fire: 0 (Minimal)

Reactivity: 0 (Minimal)

HMIS Version III Rating:

Health: 1 (Slight)

Flammability: 0 (Minimal)

Physical Hazard: 0 (Minimal)

Personal Protection: X = Consult your supervisor

HEALTH	1
FLAMMABILITY	0
PHYSICAL HAZARD	0
	X

Date of first version: 8/28/2018

### Department issuing data sheet

Contact person: see section 1: Dept. responsible for information

The information in this data sheet has been established to our best knowledge and was up-to-date at time of revision. It does not represent a guarantee for the properties of the product described in terms of the legal warranty regulations.