



# TruCal Lipid

**LOT** 34675



Cat. No. 1 3570 99 11 45 (3 x 2 mL)



**2026-04**

**In-vitro-diagnostic for veterinary use only**

Value sheet version: 1

Short Name / Analyte	Method	Reagent ID	Traceability	Value	Unit
<b>HDLC</b> HDL cholesterol	Immunoinhibition	1 3521	NIST SRM 1951b Level 2	57.5	mg/dL
				1.49	mmol/L
				0.575	g/L
	Direct, colorimetric	1 3561	Commercially available assay traceable to CDC reference method	55.6	mg/dL
				1.44	mmol/L
				0.556	g/L
<b>LDLC</b> LDL cholesterol	Select	1 4121	NIST SRM 1951b Level 2	140	mg/dL
				3.62	mmol/L
				1.40	g/L
	Direct, colorimetric	1 4131	NIST SRM 1951c Level 2	153	mg/dL
				3.96	mmol/L
				1.53	g/L
<b>Lp-PLA<sub>2</sub></b> Lipoprotein associated phospholipase A <sub>2</sub>	Colorimetric	1 7181	ε 4-nitrophenol	385	U/L
				6.41	µkat/L
<b>NEFA</b> Non esterified fatty acids	Enzymatic	1 5781	Internal Standard	56.4	mg/dL
				2.00	mmol/L
				0.564	g/L
<b>PL</b> Phospholipids	Enzymatic	1 5741	Internal Standard	287	mg/dL
				3.70	mmol/L
				2.87	g/L



DiaSys Diagnostic Systems GmbH  
Alte Strasse 9 65558 Holzheim Germany  
www.diasys-diagnostics.com